



EMERGENCY MEDICAL BAGS

Mini Bag

Art no 120224 (CPAP), 120440 (DRESSINGS), 120441 CHILD/DELIVERY)

The bags in the MiniBag series are available in three different models. The outside design and dimensions are identical while the insides of the bags differ in layout.

The MiniBag series is available in the models Medical Dressings, Child/Delivery and CPAP, and you decide which bag best suits your needs. The MiniBag is easily attached to the wall in for example an ambulance using the aluminium rail accessory.

The bag is made in the material CORDURA®, which consists of high-strength and textured polyamide 6.6 filaments, designed to withstand the harshest environments and built to last with its superior wear-resistant performance.

Accessories

Aluminium rail MiniBag, art no 112007

Care instructions

The product can be machine washed at 40°C or 60°C. Let hang dry. Do not centrifuge or tumble dry. If there is blood or vomit on the bag, rinse and wipe with bactericidal agent before washing the bag. The bag can be left open when washed.



HOW TO ORDER

Email sacci@sacci.se,
or call us at +46 (0)243-93250



PRODUCT OF SWEDEN

Weight: 900 g

Measurements (WxHxD): 320x320x120 mm

Colour: Black/neon

From Torsång to the world. With our roots in between the forests and the winding river in Torsång, the family company has delivered backpacks and carrying equipment for extreme conditions for over a hundred years. Now as then, it's the craftsmanship, innovative thinking and quality of our products that make us renowned.

Sacci Ryggsäckar AB, Gamla Tunavägen 8, S-784 36 Borlänge, SWEDEN
Phone: +46 (0)243-932 50, sacci@sacci.se

WWW.SACCI.SE

SACCI
SINCE 1914