



## EMERGENCY MEDICAL BAGS

# Respiratory care O2 Light

Art no 120670

Sacci's new respiratory care bag O2 Light holds everything you need to help a patient in need of minor or advanced respiratory assistance. The bag gives you as a user the opportunity to get a complete overview of the content quickly and easily once the bag is open.

The bag is multifunctional with an opening system that uses strong magnets to allow for quick opening and closing. The bag is fitted with ergonomically designed carrying straps and side brackets that make it possible to hang the bag on the ambulance stretcher when transporting from the accident site to the ambulance.

The bag is made in the material CORDURA®, which consists of high-strength and textured polyamide 6.6 filaments, designed to withstand the harshest environments and built to last with its superior wear-resistant performance.

**Care instructions**

The product can be machine washed at 40°C or 60°C. Let hang dry. Do not centrifuge or tumble dry. If there is blood or vomit on the bag, rinse and wipe with bactericidal agent before washing the bag. The bag can be left open when washed.

**HOW TO ORDER**

Email [sacci@sacci.se](mailto:sacci@sacci.se),  
or call us at +46 (0)243-93250



PRODUCT OF SWEDEN

**Weight:** 8,3 kg (3,2 kg without bottle)

**Measurements (WxHxD):** 330x580x250 mm

**Colour:** Black/neon

**From Torsång to the world.** With our roots in between the forests and the winding river in Torsång, the family company has delivered backpacks and carrying equipment for extreme conditions for over a hundred years. Now as then, it's the craftsmanship, innovative thinking and quality of our products that make us renowned.

Sacci Ryggsäckar AB, Gamla Tunavägen 8, S-784 36 Borlänge, SWEDEN  
Phone: +46 (0)243-932 50, [sacci@sacci.se](mailto:sacci@sacci.se)

[WWW.SACCI.SE](http://WWW.SACCI.SE)

**SACCI**  
SINCE 1914